BLOOD-BANK C

cod. 2101512

BLOOD-BANK D

cod. 2101514

BLOOD-BANK E

cod. 2101511

BLOOD-BANK F

cod. 2101513

INDICE

ESPAÑOL

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ENGLISH

12.	EQUIPMENT GENERAL INFO.....	32
12.1	Specifications	33
12.2	Packing List.....	34
12.3	Nameplate.....	34
13.	EQUIPMENT RECEPTION	34
14.	INSTALLATION.....	34
14.1	Location.....	34
14.2	Unpacking	35
14.3	Air flow	35
14.4	Levelling.....	35
14.5	Initial cleaning procedure.....	36
15.	ELECTRICAL CONNECTION.....	36
16.	STARTING UP.....	37
17.	OPERATION.....	37
17.1	Control Panel.....	37
17.2	Operation Indications.....	38
17.3	Alarms.....	38
17.4	Stirrers Connection.....	39
17.5	Electric diagram.....	40
17.6	Installing, operation & maintenance of temp. register.....	41
18.	WHEELS.....	41
19.	MAINTENANCE AND CLEANING TOPICS.....	42
19.1	Device cleaning	42
19.2	Door & hinges maintenance	43
19.3	Spare Parts & Technical Assistance	43
20.	TROUBLESHOOTING	44
21.	WARRANTY	46
22.	ANNEX (MANUAL EVJ 200).....	47

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ENGLISH

12. EQUIPMENT GENERAL INFO

This manual has been made in a simple way so that by reading it you can know the operation and maintenance of our furniture. It is recommended that you read it carefully and keep it for any reference.

The manufacturer is not responsible for any damage to people or objects that may be caused by failure to comply with the prescriptions contained in this manual. To know all the advantages of this device, please read carefully before proceeding to install it. Anyone who uses this device is recommended to read this user manual. If in doubt, consult your distributor.

This product has been manufactured under strict quality controls and meets all the requirements established by JP Selecta. Before leaving the factory, each unit has been tested to guarantee its quality. This equipment has been manufactured with recyclable materials, through an environmentally friendly production process.

This furniture complies with directive 2014/30/EC, 2014/35/EU. In addition, the CEI EN 60335-1, CEI EN 60335-2-89, EN 61000-3-2, EN 61000-6-1 and EN 61000-6-3 standards have been applied.

WARNING ! This device must be used only for the purpose described in this manual.

Electrical/electronic equipment at the end of its useful life must be managed by an authorized waste manager.

IMPORTANT: see the TECHNICAL LABEL located on the upper side of the equipment that indicates the TYPE OF REFRIGERANT.

FOR MODELS WITH R290 REFRIGERANT / HYDROCARBONS:

CAUTION – RISK OF FIRE OR EXPLOSION. FLAMMABLE REFRIGERANT. SHOULD ONLY BE SERVICED BY A QUALIFIED TECHNICIAN. DO NOT DRILL THE REFRIGERANT PIPE.

The use of electrical devices involves the implementation of basic safety instructions, such as:

- This device must be properly located and installed before installation, following the recommendations in this manual.
- Do not touch the cold surfaces of freezing devices as skin may become stuck.
- Do not store or use flammable products near the appliance.
- Unplug the appliance before any cleaning, repair or maintenance operations.

NOTE: Any manipulation of the device must be carried out by a qualified technician.

12.1 Specifications

BLOOD BANK mod. C/D are designed for fixed temperatures of +4°C (stability $\pm 1.5^\circ\text{C}$) and BLOOD BANK E/F platelet incubators are designed for fixed temperatures of +22°C. They are equipped with a hermetic compressor group mounted on an anti-vibration seat with a forced-type finned evaporator and a ventilated-type finned condenser, an air circulation turbine with a stop when the door is opened, temperature sensors housed in conditions similar to the storage bag. blood, interior LED light, double glass door, insulated and fog-proof with lock and with a circular graphic recorder.

Exterior furniture in AISI 304 stainless steel. Stainless steel interior enclosure. Reversible double glass door with easy-to-replace lock and gasket, with automatic return mechanism and manual interlock.

Interior air circulation for greater temperature homogeneity. Possibility of having an external alarm. Temp. graphic recorder. Pt100 type probe, Safety thermostat. Microprocessor temperature control that allows you to configure and view the following parameters through the display:

Security password

Set temperature

Set alarms

Calibration

Delay in starting compressor

Code	Capacity liters	Capacity bags	Internal Dimensions cm			External Dimensions cm			Consumption cooler + stirrers (W)	Refrigerant	Number of shelves	Weight (Kg)
			Height	Length	Wide	Height	Length	Wide				
2101512	560	200	153	54	66	213	69	79	497	R290 (120g)	5	132
2101514	1240	480	153	124	66	213	139	79	1034	R290 (120g)	10	190
2101511	560	200	153	54	66	213	69	79	497	R290 (120g)	5	132
2101513	1240	480	153	124	66	213	139	79	1034	R290 (120g)	10	190

12.2 Packing List

Standard equipment consists of the following components:

Description	Code	Amount
Blood-Bank C/E	2101512 / 2101511	1
Shelves	--	2
Blood-Bank D/F	2101514 / 2101513	1
Shelves	--	8
100 graphics package	0131670	1
Instructions Manual	80382	1
EC Decl. Of conformity	--	1

12.3 Nameplate

The nameplate is a label permanently affixed to the side of the equipment, which contains important electrical information, as well as data related to the refrigeration system of each unit. In addition, it includes the model and serial number.

13 EQUIPMENT RECEPTION

- All JP Selecta products are tested at the factory, evaluating their quality and performance, and not that they do not present defects.
- When you receive your appliance, it should be carefully inspected for any possible damage that may have occurred during transportation.
- If any damage is detected to the unit, you must retain all packaging material and report such damage on the carrier's delivery note. A claim must be made immediately to the transport company.
- If damage is noticed during or immediately after installation, contact your dealer immediately.

14 INSTALLATION

14.1 Location

This device is manufactured for indoor use only. Make sure that the location chosen for your equipment has adequate air circulation to ensure ventilation and efficient cooling.

Avoid locations close to heat sources, such as ovens, stoves, as well as direct solar radiation where temperatures can reach extreme values. Furthermore, a location should not be chosen in an area where the ambient humidity and temperature conditions do not exceed the specifications of the climate class of the appliance.

It must allow enough space between the equipment and the side walls, so that the 120° door opening lock can be used. Doors must be able to open a minimum of 90° in order to use the maximum available door width.

The floor of the final location must be strong enough to support the total weight of the appliance assuming it contains the maximum load of product. Additionally, it must be level and free of vibrations. Reinforce the floor if necessary.

14.2 Unpacking

The devices leave the factory on a wooden pallet and packed in resistant cardboard or wooden boxes. The box is attached to the wooden base using staples or strapping. You must first remove the staples or strapping to avoid damaging the unit when unpacking it.

All packaging materials are environmentally friendly and should be reused or recycled. Actively contribute to environmental protection by demanding recyclable packaging and nature-friendly equipment disposal methods.

NOTE: JP SELECTA does not recommend tipping the appliance forward, sideways or backwards. However, if this should occur, you must ensure that the unit remains in an upright position for at least 24 hours before connecting it, so that the compressor oil returns to the compressor.

14.3 Air flow

To ensure maximum performance of the equipment, it must be located in a location that has a continuous supply of air from both the rear and the bottom. It is recommended that the equipment be installed with respect to the wall at a separation of 75 mm of free space on each side of the unit for the same purpose. In addition, a separation from the ceiling to the furniture of 300 mm is also recommended in the case of cabinets.

A restriction in the air supply through the equipment would result in an excessive heat load on the condensing unit, which would impair its operating efficiency. At no time can the front grill of the appliance be obstructed.

NOTE: Any airflow obstruction, whether total or partial, voids the appliance warranty.

14.4 Levelling

It is very important that the appliance is perfectly level for correct operation, the doors are aligned and the unit is not subjected to undue stress.

These models are supplied from the factory with height-adjustable legs (maximum height 180 mm). In this case, you must ensure that the ground where the unit is located is level, so that it can work stably.

As an option, wheels can be supplied for all models. In case you want to install wheels, be sure to lock the front brakes so that the unit is completely stable and level. In the "Wheel Installation" section you will find detailed information about wheel adjustment.

14.5 Initial cleaning procedure

Before starting up and loading product into the appliance, you must remove the protective plastic layer that surrounds the furniture and clean it completely. If there are traces of adhesive, they must be removed with alcohol. It is recommended to clean all stainless steel surfaces of the appliance with mild soap and warm water. After cleaning, rinse with plenty of water and dry with a soft cloth.

NOTE: Never use harsh or abrasive cleaners, concentrated detergents, solvents, or chemicals to clean the equipment. Remember that products containing bleach or ammonia are very harmful to the steel surface. Avoid contamination of ferrous particles on the steel surface.

15 ELECTRICAL CONNECTION

You must check the installation voltage before connecting the equipment, checking that it is appropriate. To determine the voltage of the unit, you must check the characteristics label located on the upper exterior side of the furniture. Verify that this information exactly matches the electrical characteristics of where it will be installed.

NOTE: The device must be connected to a circuit dedicated exclusively for this purpose. Failure to comply with this requirement voids the warranty.

NOTE: The device is designed to withstand a voltage fluctuation of 5% with respect to the nominal voltage indicated on the nameplate. Compressor failure due to higher fluctuations automatically voids the warranty.

The equipment has a factory-installed hose and plug of the type shown in the following figure. If you do not have the appropriate power outlet, you must install it first.

WARNING! : If the hose or plug is damaged in any way, it could pose a serious risk. Any alteration of these components cancels the warranty.

WARNING! : JP SELECTA does not guarantee those devices connected to an extension cable (Alargos).

16 STARTING UP

Once the appliance has been installed, leveled, cleaned and electrically connected according to the instructions herein, it is ready for use. You simply have to plug it into the network.

The equipment must operate smoothly and quietly, within generally accepted standards. If you notice any unusual noise, unplug the unit immediately and check for any possible obstructions in the fans.

The device requires a certain time to reach the working temperature. You must wait for it to be reached before proceeding to load product. Continuously opening the doors hinders the equipment's ability to maintain proper cooling efficiency.

NOTE: Before loading the product, we recommend keeping the unit running for 24 hours to ensure proper operation.

NOTE: If the appliance is unplugged or disconnected, you must wait five minutes before connecting it again.

The default set temperature is 4°C (22°C in E/F models). Do not fill the equipment until the equipment reaches this temperature


17 OPERATION

17.1 Control Panel

- Large display
- Integrated keys
- Chamber probe and evaporator probe (PTC/NTC)
- Micro port input
- Alarm buzzer
- TTL MODBUS port for connection to peripherals

17.2 Operation Indications

4.1 Switching the device on and off

1. If POF = 1 (default), touch the ON/STAND-BY key for 2s.

If the device is switched on, the display will show the P5 value ("cabinet temperature" default); if the display shows an alarm code, see the section **ALARMS**.

LED	ON	OFF	FLASHING
	compressor on	compressor off	- compressor protection active - setpoint being set
	evaporator fan on	evaporator fan off	evaporator fan stop active
	cabinet light on	cabinet light off	cabinet light on by digital input
	defrost or pre-drip active	-	- defrost delay active - dripping active
	- energy saving active - low consumption active	-	-
	view temperature	-	overcooling or overheating active
	alarm active	-	-

In the annexes section, the instruction manual for the EVJ 200 programmer is attached.

17.3 Alarms

COD.	DESCRIPTION	RESET	TO CORRECT
Pr1	cabinet probe alarm	automatic	- check P0
rtc	clock alarm	manual	set date, time and day of the week
AL	low temperature alarm	automatic	check A0, A1 and A2
AH	high temperature alarm	automatic	check A4 and A5
id	open door alarm	automatic	check i0 and i1
PF	power failure alarm	manual	- touch a key - check electrical connection
C1t	compressor thermal switch alarm	automatic	check i5 and i6

The default working temperature is 4°C ±2°C in Blood Bank C/D and 22°C ±2°C in Blood Bank E/F. If the temperature goes outside the working range, whether above or below, an error occurs and is displayed on the display. See previous table. (More info in the annexes, EVJ-200)

Both the working temperature and the upper and lower limits are programmable as described in the programmer's manual in the annexes.

NOTES ON ALARMS

The acoustic signals that cause error conditions can be stopped by pressing the ON/STAND-BY key.

There is the possibility of using one of the relays on the board to activate an external alarm in case the error condition occurs (See connection diagrams).

17.4 Stirrers Connection


The equipment is factory equipped with two socket outlets with covers located in the upper right part of the cabinet, for the connection of up to 2 SELECTA shakers Mod. APL-54 / APL-108 depending on the cabinet model.

Upon request, up to 4 socket outlets can be installed to connect up to 4 APL-54 shakers.

As described in the electrical diagram, these socket outlets are protected by a magneto-thermal switch.

17.5 Electric diagram

17.6 Installing, operation & maintenance of temp. register

1. Open the register door with the key
2. Remove the needle
3. Unscrew the spindle knob and remove the paper
4. Pull the shaft to remove the recorder
5. Place the battery in the back
6. Reinsert the recorder
7. Replace the needle
8. Close the recorder door again.

Maintenance:

- Change the battery of the chart recorder every 6 months (good quality alkaline battery).
- Change the role of the recorder every week (see "Installation" section).
- Change the recorder pen every 5 years (90 meters).

18 WHEELS

These models are supplied with adjustable legs. However, there is the possibility of installing wheels optional form. To replace the legs and install the wheels, follow these steps:

Tip the appliance back carefully.

Remove the legs by unscrewing them.

Place the wheels over the same holes and screw them on.

Make sure the wheels fit properly.

Once the wheels have been checked, carefully bring the appliance to its upright position.

WARNING! Once the wheels are installed, the device must remain in a vertical position for 24 hours before being connected to ensure the return of oil to the compressor.

Screw the wheels into the same holes and make sure they are located correctly leveled

19 MAINTENANCE AND CLEANING TOPICS

19.1 Device cleaning

To clean the device, follow the following instructions:

Disconnect the equipment from the mains and remove all products inside.

Open all doors and drawers, and let the interior reach room temperature. Remove all interior accessories and clean with mild soap and warm water. Dry all accessories completely with a soft cloth.

Once the chamber has reached room temperature, clean all interior and exterior surfaces with soapy water. Rinse thoroughly and dry with a soft cloth. Not drying correctly can cause water stains to appear. Likewise, there are stainless steel cleaners that can repair and protect the protective layer of steel surfaces.

Place the accessories in their original position and connect the unit to the network.

Pitting or cracks in steel are signs of deterioration of the material. In this case, apply stainless steel cleaner capable of repairing the passivity of the steel.

NOTE: Never use steel wool, wire brushes or spatulas to clean the appliance or scouring pads.

NOTE: The cleaning products you use must be alkaline-based or chlorine-free. Any cleaner containing chlorides will damage the protective layer of the stainless steel.

Door Gasket maintenance

Door gasket requires regular cleaning to protect its elasticity, ensuring proper closure and preventing mold growth. Cleaning the gasket can be done with soapy water. Avoid using harsh cleaners and sharp utensils.


The door gasket can be easily removed by pulling it out of the door profile, and replaced by pressing it against it.

NOTE: if the gasket becomes contaminated with grease and is not cleaned, it will lose its elasticity and break in a short time.

Condenser cleaning

The condenser, located behind the rear or front grille of the unit, depending on the model, must be checked periodically; it is recommended to do so at least every 6 months. The frequency of cleaning will depend on the existing work environment. It must be ensured that air circulates freely through the condenser, so its surface must be free of dirt and grease. The dirty capacitors cause malfunctions and product loss. If the condenser coil is dirty or blocked, follow the steps detailed below.

- Disconnect the appliance from the mains.
- Remove the rear or front grille from the unit, depending on the model.
- In some models it will be necessary to remove the screws that adjust the condensing unit to the baseboard, and extract it in order to clean the condenser.
- If the condenser has a protective casing, you must unscrew it and remove it.
- Once the surface of the condenser is free, it must be cleaned using a vacuum cleaner or soft brush. Never use a metal brush.
- If the dirt is excessive, you can use compressed air to clean.
- Once clean, replace the protective casing, return the condensing unit to its original position and replace all screws.
- Finally, replace the grille and connect the appliance to the mains.

WARNING! : Never use water to clean the condenser as it could damage nearby electrical components.

19.2 Door & hinges maintenance

As time goes by and the doors are used, the hinges may move slightly. If you notice that the door is starting to misalign, you should tighten the screws that attach the hinge brackets to the cabinet.

19.3 Spare Parts & Technical Assistance

WARNING!: Make sure that the appliance is disconnected from the mains before carrying out any maintenance or repair operation. These works must be carried out by qualified personnel.

After carrying out the relevant checks, **DO NOT CARRY OUT ANY REPAIRS YOURSELF**. Contact your Technical Assistance Service, providing the model and serial number of the device (located on the characteristics label),

JP Selecta Technical Assistance Service, SAU

Tel.+34 937700877

Email:serviciotecnico@jpselecta.es

20 TROUBLESHOOTING

Many operating problems arise from causes that can be easily eliminated without the need to contact the Technical Assistance Service. The following list includes various types of problems and how to solve them.

PROBLEM	POSSIBLE REASON	POSSIBLE SOLUTION
Device not working	<ol style="list-style-type: none"> 1. The plug is not connected. 2. Electrical current is not reaching the plug because the fuse has blown or the automatic power limiter has tripped. 	<ol style="list-style-type: none"> 1. Connect the plug to the outlet current and check that electric current exists. 2. Change the fuse or reconnect the automatic power limiter.
Device is not cooling enough	<ol style="list-style-type: none"> 1. Check cutting temperatura in the controller. 2. Poorly closed door or frequent openings. 3. Obstruction of the ventilation grilles of the appliance. 4. Dirty condenser. 5. The device is exposed to direct sunlight or a heat source. 	<ol style="list-style-type: none"> 1. Reduce the cutting temperature. 2. Make sure the door does not remain open for too long. 3. Keep these areas free as indicated in the "installation" chapter of this manual. 4. Clean with compressed air or a hard hair brush (not steel). 5. Change the location of the refrigerator or protect it from these heat sources.
Noisy operation	<ol style="list-style-type: none"> 1. The appliance has not been leveled correctly. 2. Some of the inner tubes rub. 3. Loose fastening screws on some parts. 4. Fan on condenser or evaporator causing vibrations. 5. Compressor oil charge too low. 6. Loose parts in the condensing unit. 	<ol style="list-style-type: none"> 1. Level as indicated in the chapter "installation" of this manual. 2. Separate the tubes that rub. 3. Tighten the loose screws. 4. Level the equipment and tighten any loose screws. 5. If the equipment has been tipped over at any time, leave it in an upright position without connecting for 24 hours so that the oil returns to the compressor. Check for possible oil leaks.
The device is generating too much ice into the evaporator	<ol style="list-style-type: none"> 1. Poorly closed doors. 2. Excessive opening of doors. 3. Defrosting has not been carried out. 	<ol style="list-style-type: none"> 1,2. Make sure the door is not stay open for a long time. 3. Schedule a defrost.
Compressor doesnt start	<ol style="list-style-type: none"> 1. Switch open. 2. Blown fuse. 3. Defective wiring. 4. Clixon open. 5. Controller contacts open. 6. Defective relay. 7. Low gas load in the system. 	<ol style="list-style-type: none"> 1. Close the switch. 2. Replace the blown fuse. 3. Check electrical wiring. 4. Check if there is a voltage abnormally low intake current. 5. Defective controller, or device located in an area that is too cold. 6. Replace the relay. 7. Check for leaks.

Compressor starts, but stops due to overload	<ol style="list-style-type: none"> 1. Low voltage. 2. Defective unit wiring. 3. Defective starting capacitor. 4. Defective compressor. 	<ol style="list-style-type: none"> 1. Check if there is a voltage abnormally low power outlet. 2. Check wiring and electrical installation of the equipment. 3. Replace the starting capacitor. 4. Replace the compressor.
High condensing pressure	<ol style="list-style-type: none"> 1. Unit overloaded with product hot. 2. Air or non-condensable gases in the system. 3. Dirty condenser. 4. Defective condenser fan. 5. Appliance located in an area that is too hot. 6. Obstruction in expansion valve or filter. 7. Discharge valve partially closed. 8. Obstruction in discharge line. 	<ol style="list-style-type: none"> 1. Check temperature of the product, and let it cool outside the equipment if it is very hot. 2. Create a vacuum to extract air or gases. 3. Clean the condenser. 4. Replace the fan. 5. Remove equipment from very hot areas. 6. Adjust expansion valve or filter. 7. Adjust discharge valve. 8. Check and adjust discharge line.
Reduced condensing pressure	<ol style="list-style-type: none"> 1. Insufficient refrigerant charge. 2. Leaks in the system. 3. Appliance located in an area that is too cold. 	<ol style="list-style-type: none"> 1. Check if there are gas leaks refrigerant. 2. Repair system leaks. 3. Remove the equipment from very cold areas.
The compressor performs short cycles	<ol style="list-style-type: none"> 1. Differential control adjusted in intervals too small. 2. Low refrigerant charge. 3. Excessive refrigerant charge. 4. Leaking discharge valve. 	<ol style="list-style-type: none"> 1. Adjust differential with controller. 2. Check coolant pressure. 3. Adjust the coolant level. 4. Change the discharge valve. 5. Adjust refrigerant charge to avoid excess overpressures. 6. Clean the condenser.

21. WARRANTY

Dear customer, we inform you of the rules on warranty management that our company grants to its products:

- First of all, we inform you that the products manufactured and sold by JP SELECTA, SAU are capital goods intended for industrial and not domestic use. It is therefore that the guarantee applied is not regulated by the consumer and user law but by the commercial guarantee laws.
- The JP SELECTA, SAU warranty covers any manufacturing defect or any hidden defect in the device for a period of one year.
- The warranty does not cover glass breakage after delivery by JP SELECTA, SAU; nor of parts damaged by misuse or normal wear and tear.
- Any intervention on the device that affects the electrical connections, refrigeration part or electronic microcontroller not authorized by our SAT will result in the loss of the remaining warranty period on the machine.
- Any return authorized through the SAT, either for repair or replacement, is inspected at our facilities. If anomalies other than those claimed are detected that are unrelated to our manufacturing or due to misuse or wear, JP SELECTA, SAU will not be responsible for the costs of repair or replacement, which will be assumed by the customer.
- The conditions of the guarantee will not be modified unless an agreement to modify the conditions of the supply contract is previously established in writing with the client.
- **ATTENTION! NO MACHINE WILL BE ACCEPTED FOR REPAIR THAT IS NOT PROPERLY CLEANED AND DISINFECTED**


22. ANNEX (MANUAL EVJ 200)

General characteristics

- Controllers for low temperature units.
- Power supply 230 VAC.
- Incorporated clock (according to the model).
- Cabinet probe and evaporator probe (PTC/NTC).
- Door switch input.
- Compressor relay 16 A res. @ 250 VAC or 30 A res. @ 250 VAC (according to the model)
- Alarm buzzer.
- TTL MODBUS slave port for EV connect APP or BMS.
- Port for SD card data-logger module EVBD05 (according to the model).
- Models in plastic container or open-frame (according to the model)

1. User interface and main functions

1.1. Switching the device on and off

1. If POF = 1 (default), touch the ON/STAND-BY key for 2s.

If the device is switched on, the display will show the P5 value ("cabinet temperature" default); if the display shows an alarm code, see the section ALARMS.

LED	ON	OFF	FLASHING
	COMPRESSOR ON	COMPRESSOR OFF	- COMPRESSOR PROTECTION ACTIVE - SETPOINT BEING SET
	EVAPORATOR FAN ON	EVAPORATOR FAN OFF	EVAPORATOR FAN STOP ACTIVE
	CABINET LIGHT ON	CABINET LIGHT OFF	CABINET LIGHT ON BY DIGITAL INPUT
AUX 1	AUXILIARY FUNCTION 1 ON	AUXILIARY FUNCTION 1 OFF	- AUXILIARY FUNCTION 1 ON BY DIGITAL INPUT - AUXILIARY FUNCTION 1 DELAY ACTIVE
AUX 2	AUXILIARY FUNCTION 2 ON	AUXILIARY FUNCTION 2 OFF	- AUXILIARY FUNCTION 2 ON BY DIGITAL INPUT - AUXILIARY FUNCTION 2 DELAY ACTIVE

LED	ON	OFF	FLASHING
	DEFROST OR PRE-DRIP ACTIVE	-	-DEFROST DELAY ACTIVE -DRIPPING ACTIVE
	-ENERGY SAVING ACTIVE -LOW CONSUMPTION ACTIVE	*	*
	VIEW TIME	-	SET DATE, TIME AND DAY OF THE CURRENT WEEK
	VIEW TEMPERATURE	-	QUICK COOLING ACTIVE
HACCP	SAVED HACCP ALARM	-	NEW HACCP ALARMA SAVED
	ALARM ACTIVE	-	-


If Loc = 1 (default) and 30s have elapsed without the keys being pressed, the display will show the "Loc" label and the keypad will lock automatically.

1.2. Unlock keypad

Touch a key for 1s: the display will show the label "UnL".

1.3 Set the setpoint (if r3 = 0, default)

Check that the keypad isn't locked.

1.		Touch the SET key.
2.		Touch the UP or DOWN key within 15s to set the value within the limits r1 and r2 (default "-40... 50")
3.		Touch the SET key (or do not operate for 15s).

1.4 Activate manual defrost

Check that the keypad isn't locked, and that quick cooling isn't active.

1.  Touch the DEFROST key for 2s.

If P3 = 1 (default), defrost is activated provided that the evaporator temperature is lower than the d2 threshold.

1.5 Cabinet light on/off (if u1, u3, u4 or u6 = 5, default)

1.  Touch the CABINET LIGHT key.

1.6 Button-operated load on/off (if u1, u3, u4 or u6 = 9 or 10)

1.  Touch the CABINET LIGHT key (for 2s if u1, u3, u4 or u6 = 5).

If u1 or u11 = 6, the demisting switch on for the u11 duration.

1.7 Silence buzzer (if u14 = 1, default)

Touch a key.

If u1 or u11 = 11 and u9 = 1, the alarm output is deactivated.

2. Additional functions.

2.1 Activate/deactivate energy saving in manual mode

Check that the keypad isn't locked.

1.  Touch the DEFROST key.

The setpoint becomes "setpoint + r4", at maximum for HE2 duration.

2.2 View/delete compressor functioning hours

Check that the keypad isn't locked.

1.		Touch the DOWN key for 1s.								
2.		Touch the UP or DOWN key within 15s to select a label.								
	<table border="1"> <thead> <tr> <th>LAB</th> <th>DESCRIPTION</th> </tr> </thead> <tbody> <tr> <td>CH1</td> <td>view compressor functioning hundreds of hours</td> </tr> <tr> <td>CH2</td> <td>view second compressor functioning hundreds of hours</td> </tr> <tr> <td>rCH</td> <td>delete compressor and second compressor functioning hours</td> </tr> </tbody> </table>	LAB	DESCRIPTION	CH1	view compressor functioning hundreds of hours	CH2	view second compressor functioning hundreds of hours	rCH	delete compressor and second compressor functioning hours	
LAB	DESCRIPTION									
CH1	view compressor functioning hundreds of hours									
CH2	view second compressor functioning hundreds of hours									
rCH	delete compressor and second compressor functioning hours									
3.		Touch the SET key.								
4.		Touch the UP or DOWN key to set "149" (to select rCH).								
5.		Touch the SET key.								
6.		Touch the ON/STAND-BY key (or do not operate for 60s) to exit the procedure.								

2.3 View the temperature detected by the probes

Check that the keypad isn't locked.

1.		Touch the DOWN key for 1s.										
2.		Touch the UP or DOWN key within 15s to select a label.										
	<table border="1"> <thead> <tr> <th>LAB</th> <th>DESCRIPTION</th> </tr> </thead> <tbody> <tr> <td>Pb1</td> <td>cabinet temperature (if P4 = 0, 1 or 2) inlet air temperature (if P4 = 3)</td> </tr> <tr> <td>Pb2</td> <td>evaporator temperature (if P3 = 1 or 2)</td> </tr> <tr> <td>Pb3</td> <td>auxiliary temperature (if P4 = 1, 2 or 3)</td> </tr> <tr> <td>Pb4</td> <td>calculated product temperature (CPT; if P4 = 3)</td> </tr> </tbody> </table>	LAB	DESCRIPTION	Pb1	cabinet temperature (if P4 = 0, 1 or 2) inlet air temperature (if P4 = 3)	Pb2	evaporator temperature (if P3 = 1 or 2)	Pb3	auxiliary temperature (if P4 = 1, 2 or 3)	Pb4	calculated product temperature (CPT; if P4 = 3)	
LAB	DESCRIPTION											
Pb1	cabinet temperature (if P4 = 0, 1 or 2) inlet air temperature (if P4 = 3)											
Pb2	evaporator temperature (if P3 = 1 or 2)											
Pb3	auxiliary temperature (if P4 = 1, 2 or 3)											
Pb4	calculated product temperature (CPT; if P4 = 3)											
3.		Touch the SET key.										
4.		Touch the ON/STAND-BY key (or do not operate for 60s) to exit the procedure.										

3. Settings

3.1 Setting configuration parameters

1.		Touch the SET key for 4s: the display will show the label "PA".
2.		Touch the SET key.
3.		Touch the UP or DOWN key within 15s to set the PAS value (default "-19").
4.		Touch the SET key (or do not operate for 15s): the display will show the label "SP".
5.		Touch the UP or DOWN key to select a parameter.
6.		Touch the SET key.
7.		Touch the UP or DOWN key within 15s to set the value.
8.		Touch the SET key (or do not operate for 15s).
9.		Touch the SET key for 4s (or do not operate for 60s) to exit the procedure.

3.2 Set the date, time and day of the week (available in EVJ213 and EVJ214 or in EVJ203 and EVJ204 with interface EVIF25TBX connected)

	<p>N.B.</p> <ul style="list-style-type: none"> -If the device is connected to the interface EVIF25TBX, do not disconnect the device from the mains within two minutes since the setting of the time and day of the week. - If the device communicates with the APP EVconnect, the date, time and day of the week will automatically be set by the smartphone or tablet.
--	---

Check that the keypad isn't locked.

1.		Touch the DOWN key for 1s.
2.		Touch the UP or DOWN key within 15s to select the label "rtc".
3.		Touch the SET key: the display will show the label "y" followed by the last two figures of the year.
4.		Touch the UP or DOWN key within 15s to set the year.
5.	Repeat actions 3 and 4 to set the next labels.	
	LAB.	MEANING OF THE NUMBERS FOLLOWING THE LABEL
	n	month (01... 12)
	d	day (01... 31)
	h	time (00... 23)
	n	minutes (00... 59)
6.		Touch the SET key: the display will show the label for the day of the week.
7.		Touch the UP or DOWN key within 15s to set the day of the week.
	LAB.	DESCRIPTION
	Mon	Monday
	tuE	Tuesday
	UEd	Wednesday
	thu	Thursday
	Fri	Friday
	Sat	Saturday
	Sun	Sunday
8.		Touch the SET key: the device will exit the procedure.
9.		Touch the ON/STAND-BY key to exit the procedure beforehand.

3.1 Reset the factory settings

N.B.
Check that the factory settings are appropriate; see the section CONFIGURATION PARAMETERS.

1.		Touch the SET key for 4s: the display will show the label "PA".
2.		Touch the SET key.
3.		Touch the UP or DOWN key within 15s to set "149".
4.		Touch the SET key (or do not operate for 15s): the display will show the label "dEF".
5.		Touch the SET key.
6.		Touch the UP or DOWN key within 15s to set "1".
7.		Touch the SET key (or do not operate for 15s).
8.	interrupt the power supply to the device.	
9.		Touch the SET key for 2s before action 6 to exit the procedure beforehand.

4. Configuration Parameters

EVJ 203 refers to refrigeration models.

EVJ 204 refers to low temperature models.

N.	PAR.	SETPOINT	MIN... MAX.	REACH IN		TABLES					
				EVJ203	EVJ204	EVJ203	EVJ204				
	1	SP	setpoint	r1...r2	0° C	-23 °C	0° C	-23 °C			
N.	PAR.	ANALOGUE INPUTS	MIN... MAX.	REACH IN		TABLES					
				EVJ203	EVJ204	EVJ203	EVJ204				
				2	CA1	cabinet probe offset	-25... 25 °C/°F	0°C	0°C	0°C	0°C
						si P4 = 3, offset sonda aire en entrada					
3	CA2	evaporator probe offset	-25... 25 °C/°F	0°C	0°C	0°C	0°C				
4	CA3	auxiliary probe offset	-25... 25 °C/°F	0°C	0°C	0°C	0°C				

				REACH IN		TABLES	
N.	PAR.	ANALOGUE INPUTS	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204
5	P0	enable type	0 = PTC; 1 = NTC	1	1	1	1
6	P1	enable °C decimal point	0 = no; 1 = yes	0	0	0	0
7	P2	temperature unit measurement	0 = °C; 1 = °F	0	0	0	0
8	P3	evaporator probe function	0 = disabled	1	1	1	1
			1 = defrost + fan				
			2 = fan				
9	P4	configurable input function	0 = digital input	0	0	0	0
			1 = condenser probe				
			2 = critical temperature probe				
			3 = air out probe si P4 = 3, regulation temperature= product temperature (CPT)				
10	P5	value displayed	0 = regulation temperature	0	0	0	0
			1 = setpoint				
			2 = evaporator temperature				
			3 = auxiliary temperature				
			4 = air in temperature				
11	P7	inlet air weight for calculated product temperature (CPT)	0... 100 % CPT = $\left\{ \left[\frac{P7 \times (\text{inlet air T})}{100} \right] + \left[\frac{(100 - P7) \times (\text{outlet air T})}{100} \right] \right\}$	50	50	50	50
12	P8	display refresh time	0... 250 s: 10	5	15	5	15
				REACH IN		TABLES	
N.	PAR.	REGULATION	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204
13	r0	setpoint differential	1... 15 °C/°F	4°C	3°C	4°C	3°C
14	r1	minimum setpoint	-99 °C/°F... r2	-2°C	-23°C	-2°C	-23°C
15	r2	maximum setpoint	r1... 199 °C/°F	16°C	16°C	16°C	16°C
16	r3	enable setpoint block	0 = no 1 = yes	0	0	0	0
17	r4	setpoint offset in energy saving	0... 99 °C/°F	0°C	0°C	0°C	0°C
18	r5	regulation for heat or cold	0 = for cold	0	0	0	0
			1 = for heat				
19	r6	setpoint offset in overcooling	0... 99 °C/°F	0°C	0°C	0°C	0°C
20	r7	overcooling duration	0... 240 min	0 min	0 min	0 min	0 min
21	r12	position of the r0 differential	0 = asymmetric	0	0	0	0
			1 = symmetric				
				REACH IN		TABLES	
N.	PAR.	COMPRESSOR	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204
22	C0	compressor o delay afer power-on	0... 240 min	1 min	1 min	1 min	1 min
23	C1	delay between 2 compressor switch-ons	0... 240 min	1 min	1 min	1 min	1 min

N.	PAR.	COMPRESSOR	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
24	C2	compressor off minimum time	0... 240 min	1 min	1 min	1 min	1 min
25	C3	compressor on minimum time	0... 240 min	0 s	0 s	0 s	0 s
26	C4	compressor off time during cabinet probe alarm	0... 240 min	10 min	10 min	10 min	10 min
27	C5	compressor on time during cabinet probe alarm	0... 240 min	10 min	10 min	10 min	10 min
28	C6	threshold for high condensation warning	0... 199 °C/°F differential = 2°C / 4°C	70°C	70°C	70°C	70°C
29	C7	threshold for high condensation alarm	0... 199 °C/°F	80°C	80°C	80°C	80°C
30	C8	high condensation alarm delay	0... 15 min	0 min	0 min	0 min	0 min
31	C10	compressor hours for service	0... 999 h x 100 0 = disabled	0 GG	0 GG	0 GG	0 GG
32	C11	second compressor switch-on delay	0... 240 s	20 s	20 s	20 s	20 s
N.	PAR.	DEFROST	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
33	d0	automatic defrost interval	0... 99 h 0 = only manual Si d8 = 3, maximum interval	4 h	4 h	4 h	4 h
34	d1	defrost type	0... 99 h 0 = only manual Si d8 = 3, maximum interval	2	1	2	1
35	d2	threshold for defrost end	-99... 99 °C / °F	8°C	20°C	8°C	20°C
36	d3	defrost duration	0... 99 h Si P3 = 1, maximum duration	30 min	15 min	30 min	15 min
37	d4	enable defrost at power-on	0 = no; 1 = yes	0	0	0	0
38	d5	defrost delay after power-on	0... 99 min	0 min	0 min	0 min	0 min
39	d6	value displayed during defrost	0 = regulation temperature 1 = display locked 2 = dEF label	1	1	1	1
40	d7	dripping time	0... 15 min	1 min	1 min	1 min	1 min
41	d8	Defrost interval counting mode	0 = device on hours 1 = compressor on hours 2 = hours evaporator temperature < d9 3 = adaptive 4 = real time	0	0	0	0
42	d9	evaporation threshold for automatic defrost interval counting	-99... 99 °C/°F	0°C	0°C	0°C	0°C
43	d11	enable defrost timeout alarm	0 = no; 1 = yes	0	0	0	0
44	d15	compressor on consecutive time for hot gas defrost	0... 99 min	0 min	0 min	0 min	0 min
45	d16	pre-dripping time for hot gas defrost	0... 99 min	0 min	0 min	0 min	0 min

N.	PAR.	DEFROST	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
46	d18	adaptive defrost interval	0... 999 min	40 min	40 min	40 min	40 min
			if compressor on + evaporator temperature < d22				
			0 = only manual				
47	d19	threshold for adaptive defrost (relative to optimal evaporation temperature)	0... 40 °C/°F	3°C	4°C	3°C	4°C
			optimal evaporation temperature - d19				
48	d20	compressor on consecutive time for defrost	0... 999 min	180 min	180 min	180 min	180 min
			0 = disabled				
49	d21	compressor on consecutive time for defrost after power-on and overcooling	0... 500 min	200 min	200 min	200 min	200 min
			if (regulation temperature - setpoint) > 10°C/20 °F				
			0 = disabled				
50	d22	evaporation threshold for adaptive defrost interval counting (relative to optimal evaporation temperature)	-10... 10 °C/°F	-2°C	-2°C	-2°C	-2°C
			optimal evaporation temperature + d22				
51	d25	enable air out probe for defrost during evaporator probe alarm	0 = no 1 = yes	0	0	0	0
52	d26	defrost interval during evaporator probe alarm	0... 99 h	6 h	6 h	6 h	6 h
			0 = only manual				
			if d25 = 1				
				REACH IN		TABLES	
N.	PAR.	ALARMS	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204
53	A0	select value for high/low temperature alarms	0 = regulation temperature	0	0	0	0
			1 = evaporator temperature				
54	A1	threshold for low temperature alarm	-99... 99 °C/°F	5°C	10°C	5°C	10°C
55	A2	low temperature alarm type	0 = disabled	1	1	1	1
			1 = relative to setpoint				
			2 = absolute				
56	A4	threshold for high temperature alarm	-99... 99 °C/°F	10°C	10°C	10°C	10°C
57	A5	high temperature alarm type	0 = regulation temperature	1	1	1	1
			1 = evaporator temperature				
			2 = auxiliary temperature				
58	A6	high temperature alarm delay after power-on	0... 240 min	120 min	120 min	120 min	120 min
59	A7	high/low temperature alarms delay	0... 240 min	15 min	15 min	15 min	15 min
60	A8	high temperature alarm delay after defrost	0... 240 min	60 min	60 min	60 min	60 min
61	A9	high temperature alarm delay after door closing	0... 240 min	15 min	15 min	15 min	15 min
62	A11	high/low temperature alarms reset differential	1... 15 °C / °F	2°C	2°C	2°C	2°C

N.	PAR.	FANS	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
63	F0	evaporator fan mode during normal operation	0 = off 1 = on	7	7	7	7
			2 = on if compressor on				
			3 = thermoregulated (with regulation temperature + F1)				
			4 = thermoregulated (with regulation temperature + F1) if compressor on				
			5 = according to F6				
			6 = thermoregulated (with F1)				
			7 = thermoregulated (with F1) if compressor on				
64	F1	threshold for evaporator fan operation	-99... 99 °C / °F	10°C	5°C	10°C	5°C
65	F2	evaporator fan mode during defrost and dripping	0 = off 1 = on	1	0	1	0
			2 = according to F0				
66	F3	evaporator fan off maximum time	0... 15 min	10 min	10 min	10 min	10 min
67	F4	evaporator fan off time during energy saving	0... 240 s x 10	30s x10	30s x10	30s x10	30s x10
			if F0 ≠ 5				
68	F5	evaporator fan on time during energy saving	0... 240 s x 10	30s x10	30s x10	30s x10	30s x10
			if F0 ≠ 5				
69	F6	high/low humidity operation	0 = low humidity (with F17 and F18 if compressor off, on if compressor on)	0	0	0	0
			1 = high humidity (on)				
70	F7	threshold for evaporator fan on after dripping (relative to setpoint)	-99... 99 °C/°F	5°C	20°C	5°C	20°C
			setpoint + F7				
71	F8	threshold for evaporator fan operation differential	1... 15 °C/°F	2	2	2	2
72	F9	evaporator fan off delay after compressor off	0... 240 s	10 s	10 s	10 s	10 s
			if F0 = 2 or 5				
73	F10	condenser fan mode	0 = thermoregulated (with F11)	1	1	1	1
			1 = thermoregulated (with F11) if compressor off, on if compressor on				
			2 = thermoregulated (with F11) if compressor off, on if compressor on, off during defrost, pre-dripping and dripping				
74	F11	threshold for condenser fan on	0... 99 °C/°F	15°C	15°C	15°C	15°C
			differential= 2 °C/4 °F				
75	F12	condenser fan off delay after compressor off	0... 240 s	120 s	120 s	120 s	120 s
			if P4 ≠ 1				

N.	PAR.	FANS	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
76	F17	evaporator fan off time with low humidity	0... 240 s	60 s	60 s	60 s	60 s
77	F18	evaporator fan on time with low humidity	0... 240 s	10 s	10 s	10 s	10 s
N.	PAR.	DIGITAL INPUTS	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
78	i0	door switch input function	0 = disabled	5	5	0	0
			1 = compressor + evaporator fan off				
			2 = evaporator fan off				
			3 = cabinet light on				
			4 = compressor + evaporator fan off, cabinet light on				
5 = evaporator fan off + cabinet light on							
79	i1	door switch input activation	0 = with contact closed	1	1	1	1
			1 = with contact open				
80	i2	open door alarm delay	-1... 120 min	5 min	5 min	5 min	5 min
			- 1 = disabled				
81	i3	regulation inhibition maximum time with door open	-1... 120 min	15 min	15 min	15 min	15 min
			- 1 = until the closing				
82	i5	multi-purpose input function	0 = disabled	0	7 (with pressure switch in R290)	0	0
			1 = energy saving				
			2 = iA alarm				
			3 = iSd alarm				
			4 = button-operated load 1 on				
			5 = button-operated load 2 on				
			6 = device on/off				
			7 = LP alarm				
			8 = C1t alarm				
9 = C2t alarm							
83	i6	multi-purpose input activation	0 = with contact closed	1	0 (with pressure switch in R290)	1	1
			1 = with contact open				
84	i7	multi-purpose input alarm delay	0... 120 min	1	1	0	0
			if i5 = 3 or 7, compressor on delay after alarm reset				
85	i8	number of multi-purpose input activations for high pressure alarm	0... 15	0	3 (with pressure switch in R290)	0	0
			0 = disabled				
			if i5 = 3				

N.	PAR.	DIGITAL INPUTS	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
86	i9	reset counter time for high pressure alarm	1... 999 min	240 min	5 min	240 min	5 min
87	i10	door closed consecutive time for energy saving	0... 999 min	0 min	0 min	0 min	0 min
			after regulation temperature < SP 0 = disabled				
88	i13	number of door openings for defrost	0... 240	180	180	180	180
			0 = disabled				
89	i14	door open consecutive time for defrost	0... 240 min	32 min	32 min	32 min	32 min
			0 = disabled				
N.	PAR.	DIGITAL OUTPUTS	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
90	u1c	K1 output configuration	0 = compressor 1	0	0	0	0
			1 = compressor 2				
			2 = evaporator fans				
			3 = condenser fans				
			4 = defrost				
			5 = cabinet light				
			6 = demisting				
			7 = door heaters				
			8 = heater for neutral zone				
			9 = button-operated load 1 on				
			10 = button-operated load 2 on				
			11 = alarm				
12 = on/stand-by							
91	u2c	K2 output configuration	0 = compressor 1	2	4	2	4
			1 = compressor 2				
			2 = evaporator fans				
			3 = condenser fans				
			4 = defrost				
			5 = cabinet light				
			6 = demisting				
			7 = door heaters				
			8 = heater for neutral zone				
			9 = button-operated load 1 on				
			10 = button-operated load 2 on				
			11 = alarm				
12 = on/stand-by							
92	u3c	K3 output configuration	0 = compressor 1	5	5	5	5
			1 = compressor 2				
			2 = evaporator fans				
			3 = condenser fans				
			4 = defrost				
			5 = cabinet light				
			6 = demisting				
			7 = door heaters				
			8 = heater for neutral zone				
			9 = button-operated load 1 on				
			10 = button-operated load 2 on				
			11 = alarm				
12 = on/stand-by							

N.	PAR.	DIGITAL OUTPUTS	MIN... MAX.	REACH IN		TABLES		
				EVJ203	EVJ204	EVJ203	EVJ204	
✕	93	U4c	K4 output configuration	0 = compressor 1	-	2	-	2
				1 = compressor 2				
				2 = evaporator fans				
				3 = condenser fans				
				4 = defrost				
				5 = cabinet light				
				6 = demisting				
				7 = door heaters				
				8 = heater for neutral zone				
				9 = button-operated load 1 on				
				10 = button-operated load 2 on				
				11 = alarm				
12 = on/stand-by								
94	u2	enable cabinet light and button-operated load in stand-by	0 = no	0	0	0	0	
			1 = yes manual					
95	u4	enable alarm output off silencing the buzzer	0 = no	1	1	1	1	
			1 = yes					
96	u5	threshold for door heaters on	-99... 99 °C/°F	-1°C	-1°C	-1°C	-1°C	
			differential = 2 °C/4 °F					
97	u6	demisting on duration	1... 120 min	5 min	5 min	5 min	5 min	
98	u7	neutral zone threshold for heating (relative to setpoint)	-99... 99 °C/°F	-5°C	-5°C	-5°C	-5°C	
			differential = 2 °C/4 °F					
			setpoint + u7					
99	u9	enable alarm buzzer	0 = no	1	1	1	1	
			1 = yes					
				REACH IN		TABLES		
N.	PAR.	REAL TIME CLOCK	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204	
	100	Hr0	enable clock	0 = no; 1 = yes	0	0	0	0
				REACH IN		TABLES		
N.	PAR.	ENERGY SAVING	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204	
✖	101	HE2	energy saving maximum duration	0... 999 min	0'	0'	0'	0'
				- 1 = until the door opening				
				REACH IN		TABLES		
N.	PAR.	REAL TIME ENERGY SAVING	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204	
⌚	102	H01	energy saving time	0... 23 h	0 h	0 h	0 h	0 h
⌚	103	H02	energy saving maximum duration	0... 24 h	0 h	0 h	0 h	0 h
				REACH IN		TABLES		
N.	PAR.	REAL TIME DEFROST (if d8 = 4)	MIN... MAX.	EVJ203	EVJ204	EVJ203	EVJ204	
⌚	104	Hd1	1st daily defrost time	h- = disabled	-h	-h	-h	-h
	105	Hd2	2nd daily defrost time	h- = disabled	-h	-h	-h	-h
	106	Hd3	3rd daily defrost time	h- = disabled	-h	-h	-h	-h
	107	Hd4	4th daily defrost time	h- = disabled	-h	-h	-h	-h
	108	Hd5	5th daily defrost time	h- = disabled	-h	-h	-h	-h
	109	Hd6	6th daily defrost time	h- = disabled	-h	-h	-h	-h

N.	PAR.	SAFETIES	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
110	POF	enable ON/STAND-BY key	0 = no 1 = yes	1	1	1	1
111	Loc	enable keypad lock	0 = no 1 = yes	1	1	1	1
112	PAS	password	-99... 999	-	-	-	-
113	PA1	contraseña 1er Nivel	-99... 999	-	-	-	-
114	PA2	contraseña 2º Nivel	-99... 999	-	-	-	-
N.	PAR.	DATA-LOGGING EVLINK	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
115	rE0	data-logger sampling interval	0... 240 min	15 min	15 min	15 min	15 min
116	rE1	recorded temperature	0 = none, 1 = cabinet	4	4	4	4
			2 = evaporator				
			3 = auxiliary				
			4 = cabinet and evaporator				
			5 = tail				
N.	PAR.	MODBUS	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
117	LA	MODBUS address	1... 247	247	247	247	247
118	Lb	MODBUS baud rate	0 = 2,400 baud	2	2	2	2
			1 = 4,800 baud				
			2 = 9,600 baud				
			3 = 19,200 baud				
119	LP	parity	0 = none 1 = odd	2	2	2	2
			2 = even				
N.	PAR.	BLUETOOTH	MIN... MAX.	REACH IN		TABLES	
				EVJ203	EVJ204	EVJ203	EVJ204
120	bLE	enable Bluetooth	0 = no 1 = yes	0	0	0	0

5. Alarms

COD.	DESCRIPTION	RESET	TO CORRECT
Pr1	cabinet probe alarm	automatic	-check P0
Pr2	evaporator probe alarm	automatic	-check probe integrity
Pr3	auxiliary probe alarm	automatic	-check electrical connection
rtc	clock alarm	manual	set date, time and day of the week
AL	low temperature alarm	automatic	check A0, A1 and A2
AH	high temperature alarm	automatic	check A4 and A5
id	open door alarm	automatic	check i0 and i1
PF	power failure alarm	manual	-touch a key -check electrical connection
COH	high condensation warning	automatic	check C6
CSd	high condensation alarm	manual	-switch the device off and on -check C7
IA	multi-purpose input alarm	automatic	check i5 and i6
ISd	high pressure alarm	manual	-switch the device off and on -check i5, i6, i8, i9
LP	low pressure alarm	automatic	check i5 and i6
C1t	compressor thermal switch alarm	automatic	check i5 and i6
C2t	second compressor thermal switch alarm	automatic	check i5 and i6
dFd	defrost timeout alarm	manual	-touch a key -check d2, d3 and d11
FUL	SD card full alarm	manual	free up space on the SD card or replace it
Sd	No SD card inserted alarm	manual	insert the SD card or replace it