**BLOOD-BANK A                      2101519**

**BANCOS DE SANGRE  
BLOOD BANK STORAGE CABINETS**

## INDICE

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## **ENGLISH VERSION**

## GENERAL INFORMATION

- 1) *Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please notify to your J.P. Selecta supplier immediately.*
- 2) *Do not install or use the equipment without reading this handbook before.*
- 3) *This handbook must always be attached to the equipment and it must be available for all users.*
- 4) *If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.*
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) *If any modification, elimination or lack of maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.*
- 7) *Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.*

## PACKING LIST

*The standard equipment consists of the following components:*

<b>Piece</b>	<b>Part Number</b>	<b>Cant.</b>
Blood-Bank A	2101519	1
Shelves	1001864	4
Set of 100 graphs	0131670	1
Instructions manual	80403	1

## TECHNICAL FEATURES

***Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.***

Part number	Litres capac.	Bags capac.	Internal dimensions cm			External dimensions cm			Consum (W)	Power (HP)	N° shelves	Weight (Kg)
			Height	Length	Width	Height	Length	Width				
2101519	235	80	109	51	55	156	61	62	300	1/7	4	90

## EQUIPMENT DESCRIPTION

**BLOOD-BANK A** storage cabinets are designed to work at fixed temperatures at +4°C (Stability  $\pm 1,5^\circ\text{C}$ ). Equipped with an hermetic compressor mounted on shock absorbers with evaporator and forced air condenser, air cooling system which stops on opening the door, temperature sensors housed in similar conditions to the blood bag, inner LED light, double insulated anti-mist glass door and a circular chart recorder.

External cabinet in white enamelled steel and interior in painted aluminum and AISI 304 stainless steel base. Double glazed antimist, reversible lockable door, easy to change, self close door if left ajar.

Fan air circulation in order to reach homogeneous temperature distribution. External alarm ability. Circular graphic register. Pt100 probe. Safety thermostat. Temperature controlled by microprocessor which thought display allows view and control following parameters:

- Security password
- Set temperature
- Set alarms
- Calibration
- Start up cooling motor delay.

## INSTALLATION

Avoid installing the equipment near heat sources (sunny places, heating, etc). Leave a space of about 10 cm between the wall and the back of the cabinet. Place the equipment on a hard horizontal leveled floor; if it is necessary use the adjustable foot. Do not start up before 4 hours.

### VERY IMPORTANT BEFORE THE FIRST STARTING UP

Connect the equipment at mains.

The default temperature set is 4°C. Do not fill the equipment if it has not reached this temperature.

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*Be sure that the voltage supply is the same as the one indicated on the characteristics plate.*

*Do not use the equipment if it is not earthed.*

*If you have to change the plug cap, bear in mind the following:*

<b>Blue cable:</b>	<b>Neutral</b>
<b>Brown cable:</b>	<b>Phase</b>
<b>Yellow/green cable:</b>	<b>Earth</b>

### STARTING UP THE RECORDER

1. Open the recorder lid by using the key
2. Remove the needle
3. Unscrew the shaft handle and remove the paper
4. Pull out the shaft to remove the recorder
5. Insert the battery into the rear part
6. Insert the recorder again
7. Put the needle again
8. Close again the recorder lid

*Use a good quality alkaline battery to ensure a duration of more than 6 months.*

### MAINTENANCE

*Before removing the casing, disconnect the equipment from the mains.*

*The manipulation of the internal electronic circuits of the equipment by unqualified personnel can cause irreparable damage.*

*Be sure that the equipment is repaired by a J.P.SELECTA authorized technical service.*

#### BASIC MAINTENANCE

- Change recorder battery every 6 months (see "installation")
- Change recorder paper every week (see "installation")
- Change pen from recorder every 5 years (90 meters)
- Clean the evaporator of the cooling group at least once a year.

#### CLEANING

*Use alcohol with a non abrasive cloth for the cleaning of the different pieces of stainless steel and plastics.*

## MICROCOMPUTER CONTROLLER OPERATION INSTRUCTIONS

### Display and functions:

During normal operation, the controller displays the value of the probe set using parameter/4(=1 ambient default probe, =2 second probe, =3 third probe). In addition, the display has LEDs that indicate the activation of the control functions (see Table 1), while the 3 keys can be used to activate/deactivate some of the functions (see table 2).

### LEDs and associated functions

Table 1

icon	function	normal operation			start up
		ON	OFF	blink	
	Compressor	on	off	request	ON
	fan	on	off	request	ON
	defrost	on	off	request	ON
AUX	aux	output on	output off	-	ON
	alarm	all	no alarm	-	ON
	clock	RTC fitted and enabled, at least 1 time band set	RTC not fitted or disabled, not even 1 time band set	-	ON if RTC fitted

Table of functions activated by the buttons-models S, X, Y & C

Table 2

button		normal operation		start up	
		pressing the button alone	pressed together		
	up ON/OFF	more than 3s toggle ON/OFF	pressed together start/stop continuous cycle	-	
	down defrost	more than 3s: start/stop defrost		Pressed together start parameter reset procedure	for 1s display firmware vers.code
	set mute	-1s: display/set the set point -more than 3s: access parameter Setting menu(enter password 22) -mute audible alarm (buzzer)			for 1s RESET current EZY set

**Setting the set point (desired temperature)**

- Press SET for 1 sec the set value will start flashing after a few moments.
- Increase or decrease the value using UP or DOWN.
- Press SET to confirm the new value.

**Switching the device ON/OFF**

Press UP for more than 3 sec. The control and defrost algorithms are now disabled and the instrument displays the message "OFF" alternating with the temperature read by the set probe.

**Manual defrost (models S, X, Y and C only)**

Press for DOWN more than 3 sec (the defrost starts only the temperature conditions are valid).

**Continuous cycle (models S, X, Y and C only)**

Press UP and DOWN together for more than 3 sec.

**Access and setting type F (frequent) and type C (configuration) parameters**

1. Press SET for 3 sec (the display will show "PS")
2. To access the type F and C parameter menu, enter the password "22" using UP/DOWN;  
To access the F parameter menu only, press SET (without entering the password);
3. Scroll inside the parameter menu using UP/DOWN;
4. To display/set the values of the parameter displayed, press SET, then UP/DOWN and finally SET to confirm the changed (returning to the parameter menu).  
To save all the new values and exit the parameter menu, press SET for 3 sec.  
To exit the menu without saving the changed values (exit by timeout) do not press any button for at least 60 sec

Table of alarms

Alam code	buzzer and alarm relay	LED	Description	Parameters involved
E0	active	ON	probe 1 error=control	-
E1	inactive	ON	probe 2 error=defrost	[d0=0/1]
E2	inactive	ON	probe 3 error=condenser	[A4=10]
IA	active	ON	external alarm	[A4=10][+A7]
dOr	active	ON	open door alarm	[A4=7/8][+A7]
LO	active	ON	low temperature alarm	[AL][Ad]
HI	active	ON	high temperature alarm	[Ah][Ad]
EE	inactive	ON	unit parameter error	-
EF	inactive	ON	operating parameter error	-
Ed	inactive	ON	defrost ended by timeout	[dP][dt][d4][A8]
dF	inactive	ON	defrost running	[d6=0]
cht	inactive	ON	condenser dirty pre-alarm	[A4=10]
CHT	active	ON	condenser dirty alarm	[A4=10]
EtC	inactive	ON	clock alarm	

Our products have been precisely modified before leaving factory. To avoid damaging compressor unit or other malfunctions, users must not modify the microcomputer parameters privately.

**NOTES ON ALARMS**

- The acoustic signals that cause the error conditions can be stopped by pressing the ON/STAND-BY key.
- There is the possibility of using one of the board relays to activate an external alarm in the event that the error condition occurs (See connection diagrams).

## DESVIATION OF TEMPERATURE

*The default working temperature is 4°C ±2°C in Blood Bank A. If the temperature rises outside the working range, either above or below, an error occurs which is displayed on the screen (See previous table).*

*Both the working temperature and the upper and lower limits are programmable as explained in the section "Modification of parameters".*

## GUARANTEE

*This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.*

*Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.*

### Notice to customers:

*The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.*