PLAC-CENTER

4000501

**PLACA REFRIGERANTE
COOLING PLATE**

INDICE

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GENERAL INFORMATION

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY EQUIPMENT TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 2009/104/CE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the equipment with liquids which can give off vapours capable of making explosive mixtures.

PACKING LIST

The standard equipment consist of the following components:

Description	Code	Quantity
Plac-center	4000501	1
Instruction Manual	80052	1

TECHNICAL FEATURES

Voltage supply 115-230V 50/60 Hz according to the characteristics plate indications.

Code		4000501
Plate dimensions cm.	Length	30
	Width	15
Overall dimensions cm	Height	27
	Length	33
	Width	50
Compressor power HP		1 / 3
Consumption W		160
Gas type		R 134A
Gas load gr.		150
Weight Kg.		15

EQUIPMENT DESCRIPTION

The cooling plate **PLAC-CENTER** is specially suitable for maintaining substances in solid state.

For adjustable temperatures on plate from 5°C up to -10°C, with electronic regulation and digital reading.

AISI 304 stainless steel outer casing and cooling plate in black anodized aluminium. Hermetically-sealed compressor installed on an anti-vibration mounting with ventilated evaporator, regulating valve and condenser.

Control panel equipped with main switch with indicator lamp, analogue electronic temperature regulator and digital reading (2 digits) with PTC probe and resolution ± 1 digit.

INSTALLATION

Place the cooling plate on a flat horizontal, level surface, leaving a space at the back and sides of the apparatus of about 10cm.

CAUTION!!! IMPORTANT FOR YOUR SAFETY!

Do NOT use the plate for thermic treatment of products which made form explosive gases, not even in special protection aerea.

Be sure that the voltage supply is the same as the one indicated on the characteristics plate.

Do not use the apparatus if it is NOT earthed.

If you have to change the plug bear in mind the following:

Blue cable: Neutral.


Brown cable: Phase.

Yellow/green cable: Earth.

OPERATION

CONTROL PANEL:

1. Main switch with indicator lamp.
2. Compressor operation indicator lamp.
3. Temperature controller.

STARTING UP:

Before the introduction of the first load, it is advisable to leave the equipment working for about two hours at the preset temperature.

OPERATION:

1. Connect the plate to the mains.
2. Press the main switch (1).
3. Adjust set point (SP).

SET key (5): press for 5 seconds to vary the SP (Set Point). Modify with the **Up** (6) and **Down** (7) keys to select the setpoint temperature. Validate by pressing the **SET** key (5) again. The stabilization time at the desired temperature will vary depending on the ambient temperature.

Up key (6): press for 5 seconds to start or stop defrosting. In the programming menu it is possible to move through the different levels, and to vary a parameter when adjusting it.

Down key (7): press for 5 seconds to activate Stand-by mode. To return to normal mode, press again for 2 seconds. In Stand-by mode, the equipment does not perform any action and only the indicator (1) is shown on the screen.

In the programming menu, it allows you to scroll through the different levels, and to vary a parameter when adjusting it.

Indicator leds and keys:

1. Stand-by mode indicator
2. Programming mode indicator
3. Display
4. Active Cooling Indicator
5. SET key
6. Up / Defrost key
7. Down / Stand-by key

MAINTENANCE

Before removing the side covers, disconnect the equipment from the mains.

The guarantee will not have any effect if the equipment is manipulated by somebody who has not been authorized by J.P.SELECTA, s.a.u.

CLEANING:

Remove the metallic casing and clean the condenser wings with the help of a brusher or compressed air. It is recommended to do it periodically.

GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by an incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by unauthorized personnel by J.P. SELECTA, s.a.u. cancels the guarantee automatically.

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.